
BRAUN

Braun No Touch+ touch thermometer
百靈免接觸及接觸額溫槍

No touch + forehead mode
免接觸式及接觸額式

AgeSmart age-adjustable fever guidance

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Also measures food, object or liquid temperature
可測量食物，物體或液體溫度


How does it work 如何運作

The No Touch+ touch thermometer makes temperature taking easier, faster and stress-free.

The unique dual technology means that you can take temperature two ways - by placing the thermometer on the forehead, and also by holding the thermometer in front of the forehead using the no touch mode. Designed for the whole family, with no touch mode being ideal to conveniently monitor temperature without having to disturb your child and touch mode for comfortable self-measurement.

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
The age-adjustable fever guidance that takes the guesswork out of interpreting temperature for the whole family.

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
Clinical research shows that the definition of fever changes with a child's age*. Many parents and guardians are not aware of this. To help them interpret a child's temperature at first sight, Braun has unique AgeSmart technology.

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Simply select your child's age and the color-coded display will help you understand the results and act accordingly - green (normal temperature), yellow (moderate fever) and red (high fever).

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How to use? 如何使用?

Taking an individuals temperature 測量體溫

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
4  Position between eyebrows and hold on forehead OR up to 2.5 cm (1 inch) away from forehead.

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 Yellow light will show you where you are aiming.

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Taking a food & bath temperature 測量食物或洗澡水溫度

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3  Press the food and bath button and water drop icon appears on display screen.

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
4  Position up to 6.35 cm (2.5 inches) from food or bath water. Do not submerge.

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NOTE: When measuring food or bath water temperatures, the thermometer must be aimed directly at the food or bath water. Keep out of direct line of the steam, steam will impact the reading by fogging over the sensor. Do not aim at the sides of the container as this will impact the measurement.

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5  Press temperature button
按下測量溫度按鈕

NOTE: Each time you power on, the thermometer will default to human temperature mode.

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What is fever? 何謂發燒?

If your body temperature is consistently above your normal temperature range*, you are said to have fever. Fever is a defense mechanism the body uses to fight infections. However, fever alone isn't an illness, it's a symptom.

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* Practice taking temperatures on yourself to determine the normal temperature range.

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Temperature taking hints 測量體溫小秘訣

- Let the thermometer and individual being measured acclimate to the temperature of the room in which the measurement will be taken.
- Thermometer needs to acclimate to the room temperature for 10 minutes before taking a measurement.
- The individual being measured should be indoors for 30 minutes before taking a measurement. External temperature can impact skin temperature.
- 讓額溫槍及接受測量體溫的人士先適應擬進行測量的房間的室溫。
- 免接觸額溫槍須在室內放置 10 分鐘才可操作。
- 要測量體溫的人士須在室內逗留 30 分鐘才可測量體溫。外部環境溫度會影響皮膚溫度。

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

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© 2020 All rights reserved.