BRAUN

ThermoScan® 7  
Infrared ear thermometer  
百靈耳溫槍

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BRAUN

Do you know the  
definition of fever  
changes with age?  
My thermometer does.

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Clinical research shows that the definition  
of fever changes with a child's age<sup>3</sup>. Many  
parents and guardians are not aware of  
this. To help them interpret a child's  
temperature at first sight, Braun has  
developed unique AgeSmart™ technology.

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Braun AgeSmart™ technology provides  
age-adjustable fever guidance to take  
the guess work out of interpreting  
temperature for the whole family.  
百靈 AgeSmart™ 採用可隨年齡調節的發燒  
指示技術，讓您為家人測溫時毋須再猜測是  
否發燒。


Fever guidance  
發燒指示功能

Large color-coded display provides  
instant fever detection for the whole  
family.  
顯示幕以不同顏色顯示體溫狀況，讓您輕易  
解讀各家庭成員的體溫。

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Backlit display  
背光功能

Helps take a temperature even  
when it's dark.  
能在黑暗環境下測量及讀取體溫。

Mx9

Memory function  
記憶功能

Recalls last 9 temperature readings.  
可保存最近 9 次讀數。

Safe and hygienic  
使用安全及衛生

Braun disposable lens filters help to  
eliminate the risk of cross-contamination  
and ensure accurate readings.  
百靈一次性耳套可避免交叉感染和確保測溫的  
準確度。

O+

Suitable for newborns  
適合初生嬰兒

Pre-warmed tip  
測溫準確

Pre-warmed tip technology is designed to  
minimize the cooling effect of the tip inside  
the ear canal, providing accurate  
measurements time after time.  
預熱測溫頭可將測溫頭在耳道內產生的冷卻效應降  
低，確保連續測溫的準確度。

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
3. Herzog L, Phillips S. Addressing Concerns About Fever. Clinical Pediatrics.2011; 50(#5):  
383-390.  
4. Nimah, Marianne M., et. al. "Infrared tympanic thermometry in comparison with other  
temperature measurement techniques in febrile children." Pediatr Crit Care Med. 7:1 (2006):  
48-55.


Braun's Ear Thermometer  
Braun 百靈耳溫槍

Braun ThermoScan® measures the infrared heat  
generated by the eardrum and surrounding tissue.  
Because the eardrum shares blood supply with the  
temperature control centre in the brain, measurement  
with Braun ThermoScan® accurately reflects changes  
in core body temperature.  
ExacTemp™ stability indicator uses a light and beep to  
confirm proper positioning. Pre-warmed tip  
technology is designed to minimize the cooling effect  
of the tip inside the ear canal, providing accurate  
measurements time after time.

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Certain trademarks used under license from The Procter & Gamble Company or its affiliates. ThermoScan® and ExacTemp™ are trademarks of Helen of Troy Limited and/or its affiliates. 部份商標為寶潔公司或其關聯公司認可使用。 ThermoScan® 和 ExacTemp™ 為 Helen of Troy Limited 及/或其關聯公司的註冊商標。 © 2025, All rights reserved.

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



Age-adjustable fever guidance for the whole family.  
隨著年齡組別而調節的發燒指示功能。

Always replace disposable lens filters to maintain accuracy and hygiene.

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### What is fever? 何謂發燒?

If your body temperature is consistently above your normal temperature range\*, you are said to have fever. Fever is a defense mechanism the body uses to fight infections. However, fever alone isn't an illness, it's a symptom. 當身體體溫高於正常的體溫範圍\*，便稱為發燒。發燒是身體為了打敗外來感染而產生的防護機制。發燒並不是一種疾病，而是其他疾病的徵兆。

### Would the tip on the thermometer affect accuracy reading? 測溫頭會否影響體溫的準確度?

Yes. The tip on other ear thermometers can cool the area where the reading is taken, leading to a higher risk of inaccurate readings. The Braun pre-warmed tip is designed to minimize the cooling effect of the tip inside the ear canal, providing accurate measurements time after time. 會。市面上耳溫槍產品的測溫頭，在測量時，會令測溫部位溫度降低，從而可能導致測溫結果失準。百靈預熱測溫頭，在使用前已做預熱處理，可將測溫頭在耳道內產生的冷卻效應降低，有助確保連續測溫的準確度。

### How can parents reduce the risk of infection when their children have a fever? 小孩發燒時，父母有甚麼方法能減低交叉感染的機會?

Infants and children can suffer from fever many times throughout the year, and parents usually take their child's temperature several times during each bout of fever, depending on how unwell the child feels. Taking measurement with Braun disposable lens filter enable hygienic use without having to clean or disinfect the sensor and virtually eliminate cross-contamination. What's more, fitting a new Braun lens filter each time ensures that dirt or ear wax does not build up and clog the filter membrane, which could lead to inaccurate temperature readings.

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